










HUGE Digitalife™ is our state-of-the-art Digital Denture Design Concept. It features the harmonious integration of prefabricated teeth and digitally generated denture bases, alongside permanent zirconia restorations. Digitalife™ Solution ensures precision, durability, and superior aesthetics in creating complete, functional, and aesthetically pleasing dental solutions.

## >> DIGITAL PRINT DENTURE




 <b>01 &gt;&gt; Oral Scanning</b> i-Vinci Intraoral Scanners	 <b>02 &gt;&gt; CAD Process</b> HUGE Tooth Libraries available in <i>sshape</i> and <i>exocad</i>	 <b>03 &gt;&gt; 3D Printing</b> Print the denture base
 <b>06 &lt;&lt; Satisfactory Result</b> Bring Convenient Smiles	 <b>05 &lt;&lt; Digital Bonding</b> HUGE Self Curing Resin	 <b>04 &lt;&lt; Teeth Preparation</b> HUGE Synthetic Polymer Teeth

Impression	i-Vinci Intraoral Scanner PERFIT Silicone Impression Materials	
Design	Digitalife™ Dental Cloud (International)	
Printed Denture Base	NOBILCAM® Digitalife™ 3D Print Denture Base Resin NOBIL-3D Printer/ NOBILCAM® 3D Printer	
Denture Teeth	NOBILCAM® Digitalife™ 3D Print Prime Tooth Resin(Planning) NOBILCAM® Digitalife™ 3D Print Selected Tooth & Trial Placement Resin(Planning) Sonning Digitalife™ Version Denture Teeth MAIST Digitalife™ Version Denture Teeth(Planning) KAILI Digitalife™ Version Denture Teeth (Planning)	

## >> DIGITAL MILL DENTURE




 <b>01 &gt;&gt; Oral Scanning</b> i-Vinci Intraoral Scanners	 <b>02 &gt;&gt; CAD Process</b> HUGE Tooth Libraries available in <i>sshape</i> and <i>exocad</i>	 <b>03 &gt;&gt; Milling</b> Milling the denture base
 <b>06 &lt;&lt; Satisfactory Result</b> Bring Convenient Smiles	 <b>05 &lt;&lt; Cementing</b> Bonding the denture teeth in the denture base	 <b>04 &lt;&lt; Teeth Preparation</b> Milling the dental arch or using prefabricated teeth

Impression	i-Vinci Intraoral Scanner PERFIT Silicone Impression Materials	
Design	Digitalife™ Dental Cloud (International)	
Milled Denture Base	Brilliance PMMA Discs NOBILCAM®-5X Milling Machine	
Denture Teeth	Zenith & Radiance PMMA Disc Sonning Digitalife™ Version Denture Teeth MAIST Digitalife™ Version Denture Teeth (Planning) KAILI Digitalife™ Version Denture Teeth (Planning)	

## >> ZIRCONIA RESTORATION (CROWN, BRIDGE AND FULL-ARCH)



 <b>01 &gt;&gt;</b> Precise impression-taking	 <b>02 &gt;&gt;</b> 3D Print Dental Model	 <b>03 &gt;&gt;</b> Milling the denture arch
 <b>06 &lt;&lt;</b> Bonding System	 <b>05 &lt;&lt;</b> Re-sinter and Polish	 <b>04 &lt;&lt;</b> Sinter and Glaze

Impression	i-Vinci Intraoral Scanner PERFIT Silicone Impression Materials	
Design	Digitalife™ Dental Cloud (International)	
Temporary restoration	Zenith & Radiance PMMA Disc	
Milled zirconia	Zirking Zirconia Disc (MM4D specially for full-arch) NOBILCAM®-5X Milling Machine	
Bonding	HugeBond Universal FlipPro Light Cure Dental Adhesive TopCEM Dual Cure Resin Cement	

