



HUGE Your Partner
in Prosthetic
Dentistry






Digitalife™ Tooth Mould Chart Denture Solutions

Designed for Digitalife™ Denture Solution



Shandong Hugu Dental Material Corporation
Add / No.68 Shanhai Road, Donggang District, Rizhao City, Shandong Province,
276800, P.R. China.
Tel / +86 (633) 2277268 marketing@hugedental.com www.hugedental.com

-  Huge Dental
-  Huge Dental
-  Huge Dental

Step Ahead:

Digitalife™ teeth are specifically designed to fit with digitally manufactured denture base, offering a refined solution for HUGE traditional synthetic polymer teeth. After parameter and moulds adjustments, Digitalife™ teeth integrate more effectively with digital denture bases in 70% of cases.

Sonning is proud to be the first carded teeth to offer a Digitalife™ version.



Designed for Digitally Fabricated Dentures



Advanced Mould Design

Sonning Digitalife contains selected 11 anterior and 6 posterior moulds that are designed for the digital denture workflow. They can meet up to 70% of cases needs thanks to its shorter neck, thinner back and undercut-free basal design. Elevating your digital denture workflow!



Favored Esthetics

Inspired by *Sonning* morphology which has gained much affection from users, Sonning Digitalife Version balances esthetics with functionality, providing a more vivid look, while having a lovely form, making it a good choice for a wider range of patients.



Elevated Efficiency


Sonning Digitalife is compatible with our Digitalife™ Denture Solutions, which saves on labor and material costs while enabling rapid bonding with the denture base.



Deep Integration with Digitalife™ Denture Solution

Sonning Digitalife has pre-set libraries include combinations for Physiological, Gerber, Lingualized and Crossbite for all 3Shape Dental System users.

Well-designed Tooth Moulds



Elevating your Denture Design

With four preset occlusal concepts in denture tooth library, Sonning Digitalife can meet up to 70% cases' needs and saves on your labor cost.



Crossbite

Gerber

Lingualized

Physiological

HUGE Digitalife™ Denture solution

HUGE Digitalife™ is our state-of-the-art Digital Denture Design Concept. It features the harmonious integration of base and denture teeth.

We've now made our entire range of tooth libraries available for 3Shape and exocad. With our functional designs and aesthetically pleasing results, we hope both your professional team and your patients will appreciate the quality of dentures made from digitally fabricated bases and prefabricated teeth.


What can you expect:

- Time-saving process* ✓
- Aesthetic appearance* ✓
- Precise result* ✓
- Meet your individual needs* ✓
- Deep integration with digital dentistry* ✓




01 >> Oral Scan

i-Vinci Intraoral Scanners




02 >> CAD Process

HUGE Tooth Libraries available in 3shape and exocad



03 >> Print Base

3D Denture Base Resin 3D Printer



05 << Great Result

Enjoy Perfect Smiles



04 << Teeth Bonding

*HUGE Prefabricated Synthetic Polymer Teeth
HUGE Digital Bond Kit*



03 << Mill Base

HUGE Pink PMMA BLOCK Milling Machine

  **SQUARE**



  **OVOID**




  **TAPERING**




LOWER ANTERIORS



S4 All measurements in millimeters.



POSTERIORS 28°



COMBINATIONED ARTICULATIONS

Upper anteriors	Lower anteriors	Posterior 28°
S4	L7	32
O5	L8	34
T3	L6	32
T5	L7	32
T6	L7	32
T9	L11	32

Sonning
Digitalife Version Adapts to Digitalife™ Denture Solutions


Focus on aesthetics, function and performance.

Moulds

- >> 6 upper anteriors moulds
- >> 5 lower anterior moulds
- >> 6 posterior moulds

Shades


- >> Classical A2 & A3



Packaging

- >> Anterior: 6×1×16/Box (6pcs/card, 16 cards/box)
- >> Posterior: 8×1×12/Box (8pcs/card, 12 cards/box)

(Tooth Moulds Digitally Available)



3shape 
(available now)